



Email : sales@laboratoryinstrumentindia.com

Product Name :
Stirred Tanks In Series

Product Code :
LBNY-0005-12700014



Description :

Chambered bottom made of stainless steel as heat exchanger for connection
3 identical stirred tank reactors made of glass connected in series
Sensor for measuring the conductivity and temperature
Stirred tanks in series for connection to supply unit
Display of conductivity and temperature via
Delivery between stirred tanks via 2 peristaltic pumps of the supply unit
Temperature control in the reactor via
Small reactor capacity for less consumption of chemicals.

Technical Specification :

3 Reactors
Inside diameter: each 70mm
Outer diameter: each 80mm
Height: each 140mm
Speed stirrer: 3x approx. 330min⁻¹
Reactor capacity: each approx. 350ml
Measuring Ranges:
Temperature: 0...60°C
Conductivity: 0...100ms/cm.



Laboratory Instrument India