



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Stirred Tanks In Series

**Product Code :**  
LBNY-0005-12700014



#### **Description :**

Chambered bottom made of stainless steel as heat exchanger for connection  
3 identical stirred tank reactors made of glass connected in series  
Sensor for measuring the conductivity and temperature  
Stirred tanks in series for connection to supply unit  
Display of conductivity and temperature via  
Delivery between stirred tanks via 2 peristaltic pumps of the supply unit  
Temperature control in the reactor via  
Small reactor capacity for less consumption of chemicals.

#### **Technical Specification :**

3 Reactors  
Inside diameter: each 70mm  
Outer diameter: each 80mm  
Height: each 140mm  
Speed stirrer: 3x approx. 330min-1  
Reactor capacity: each approx. 350ml  
Measuring Ranges:  
Temperature: 0...60°C  
Conductivity: 0...100ms/cm.



**Laboratory Instrument India**