



Email : sales@laboratoryinstrumentindia.com

Product Name :
Plug Flow Reactor

Product Code :
LBNY-0005-12700015



Description :

Straight Glass Tube With Fixed Bed From Glass Spheres As Reactor
Transparent Double Jacket From Pmma For Cooling And Heating
Plug-Flow Reactor For Connection To Supply Unit
Sensor For Measuring The Conductivity And Temperature Via
Air Vessel For Damping Of Pulsation
T-Piece With Nozzle For Mixing The Reactants
Temperature control in the reactor via.

Technical Specification :

Water Bath
Material: Pmma
Capacity: Approx. 0,4l
Inside Diameter: 70mm.
Plug Flow Reactor
Material: glass
Inside Diameter: 40mm
Height: 400mm.

Laboratory Instrument India

