



Email : sales@laboratoryinstrumentindia.com

Product Name :
Continuous Stirred Tank Reactor

Product Code :
LBNY-0005-12700016



Description :

Chambered bottom made of stainless steel as heat exchanger for connection to
Sensor for measuring the conductivity and temperature via
Height-adjustable overflow for changing the reactor volume
Temperature control in the reactor via
Continuous stirred tank reactor for connection to supply unit
Glass tank
Reactor with stirrer.

Technical Specification :

Stirred Tank Reactor
Height: 120mm
Outer diameter: 110mm
Inside diameter: 100mm
Speed stirrer: approx. 330min⁻¹
Adjustable volume: 270...750ml.



Laboratory Instrument India

