



Email : sales@laboratoryinstrumentindia.com

Product Name :
Ion Exchange

Product Code :
LBNY-0005-12700021



Description :

Tank With 4 Chambers For Raw Water, Rinsing Water, Acid And Caustic
Diaphragm Pump To Transport Raw Water, Rinsing Water, Acid And Caustic
Softening And Desalination With Ion Exchange
Collecting Tank For Treated Water, Rinsing Water, Acid And Caustic
Continuous measurement of conductivity and flow rate
Cation And Anion Exchangers Usable Separately And In Combination
Regeneration Of Ion Exchangers.

Technical Specification :

Ion Exchange



Laboratory Instrument India