




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Microcontroller Trainer

**Product Code :**  
LBNY-0002-610005



**Description :**

Microcontroller Trainer

**Technical Specification :**

A full-featured development board for high pin-count 8-bit Microchip PIC microcontrollers. It contains many on-board modules necessary for the development variety of applications, including graphics, Ethernet, USB, data logging and others. Designed for use with surface-mount package versions of the PIC MCU, mounted onto an MCU card. Arduino/Click Boards

Microcontroller Trainer, Microcontroller Trainer Equipments, Microcontroller Trainer Tools, Microcontroller Trainer Tool Kits, Microcontroller Trainer Manufacturers, Microcontroller Trainer Suppliers from India, China, Kenya



**Laboratory Instrument India**