



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Track String Adapter

**Product Code :**  
LBNY-0004-780021



**Description :**

Track String Adapter

**Technical Specification :**

Mechanical bench accessory 19:

Track string adapter:

The Dynamics Track Mount is used to mount the Rotary Motion Sensor to the Dynamics Track, allowing it to act as a high resolution, bi-directional Smart Pulley.

When used in conjunction with the Dynamics Track Mount, the Track String Adapter allows a Rotary Motion Sensor to continuously monitor the Dynamics Cart position.

To attach on the rotary motion sensor pulley and the cart via a clip.

Include:

A bracket with pulley

A cart string clip

A thread



**Laboratory Instrument India**

