



Email : sales@laboratoryinstrumentindia.com

Product Name :
Adsorption

Product Code :
LBNY-0005-12700023



Description :

Safety adsorber for closed water circuit
Continuous process
2 adsorbers with activated carbon filling
Water temperature control
Digital temperature indication
Adsorber with 8 sampling points
Metering pump for concentrated adsorbate solution
Change of adsorbate concentration and contact time
Pump for recirculating the treated water
Flow rate adjustable.

Technical Specification :

Adsorber and safety adsorber
Capacity: each 1700cm³
Inside diameter: each 60mm
Height: each 600mm
Circulation pump
Max. Flow rate: 180l/h
Max. head: 10m
Tanks
Adsorbate solution: 45L
Treated water: 45L

Heater
Max. power: 500W
Metering pump
Max. Flow rate: 2,1l/h
Max. head: 160m.



Laboratory Instrument India