



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Adsorption

**Product Code :**  
LBNY-0005-12700023



#### **Description :**

Safety adsorber for closed water circuit  
Continuous process  
2 adsorbers with activated carbon filling  
Water temperature control  
Digital temperature indication  
Adsorber with 8 sampling points  
Metering pump for concentrated adsorbate solution  
Change of adsorbate concentration and contact time  
Pump for recirculating the treated water  
Flow rate adjustable.

#### **Technical Specification :**

Adsorber and safety adsorber  
Capacity: each 1700cm<sup>3</sup>  
Inside diameter: each 60mm  
Height: each 600mm  
Circulation pump  
Max. Flow rate: 180l/h  
Max. head: 10m  
Tanks  
Adsorbate solution: 45L  
Treated water: 45L

---

Heater  
Max. power: 500W  
Metering pump  
Max. Flow rate: 2,1l/h  
Max. head: 160m.



**Laboratory Instrument India**