



Email : sales@laboratoryinstrumentindia.com

Product Name :
Advanced Oxidation

Product Code :
LBNY-0005-12700024



Description :

Flow rate adjustable
Measurement of temperature and flow rate
Digital temperature indication
Advanced oxidation process
Falling film reactor with uv lamp
Discontinuous operation
Use of hydrogen peroxide and uv light
Formation of hydroxyl radicals (oh radicals)
Protection device against UV radiation.

Technical Specification :

UV lamp
Power: 120W
Emitted wavelength: 254nm
Falling film reactor (tube)
Material: glass
Diameter: 130mm
Height: 1000mm
Tank
Capacity: 10L
Pump
Max. head: 9m

Max. Flow rate: 360l/h.



Laboratory Instrument India