



Email : sales@laboratoryinstrumentindia.com

Product Name :
Convection Drying

Product Code :
LBNY-0005-12700025



Description :

Air Heating With Controlled Heater
Digital Balance For Measuring The Change Of Weight During Drying
Drier For Investigating Convection Drying Of Solids
Adjustment Of Air Velocity Via Speed Of Fan
1 Combined Sensor For Measurement Of Humidity And Temperature Before And After The Solid Sample Drying On 4 Corrosion Resistant Plates In A Drying Channel With An Air Flow
Digital stopwatch, battery operated
1 Air Velocity Sensor.

Technical Specification :

Drying channel
4 drying plates: 398x320mm
Length: 2340mm (with fan)
Internal dimensions: 350x350mm
Heater
With adjustable temperature limiter
Power: 0...6750w
Fan
Power: 33W
Max. speed: 950min-1
Max. Output: 700m/h
Balance

Application temperature: 075C
Measuring range: 010000g
Resolution: 0,1g.



Laboratory Instrument India