



Email : sales@laboratoryinstrumentindia.com

Product Name :
Discontinuous Rectification

Product Code :
LBNY-0005-12700031



Description :

Packed column with raschig rings
Vacuum mode possible with water jet pump
Discontinuous rectification with packed and sieve tray column
All tanks made of duran glass and stainless steel
Adjustment of reflux ratio using valves
Electrically heated evaporator
Tank for top product
Interchangeable columns
Sieve tray column with 8 trays
Condenser and phase separation tank for top product.

Technical Specification :

Water jet pump: final vacuum: approx. 200mbar
Columns: internal diameter: 50mm, height: 765mm
Evaporator
Power output: 0...4kw
Tank: approx. 10L
Tanks
Top product: approx. 2000ml
Phase separation: approx. 500mL
Heat transfer surface
Top product condenser: approx. 0,04m².



Laboratory Instrument India