



Email : sales@laboratoryinstrumentindia.com

Product Name :
Suspension Production Unit

Product Code :
LBNY-0005-12700034



Description :

Eccentric Screw Pump, With Pressure Cut-Out Switch, Dry-Running Protection And Adjustable Speed, To Deliver the suspension
Supply Unit To Produce And Deliver Suspensions For Experimental Filtration Units
Stirred Tank With Lid And Stirring Machine To Prepare A Suspension

Technical Specification :

Tank: 200L, stainless steel
Pump
Manometer measuring range: 0...10bar
Max. Head: 50m
Max. Flow Rate: Approx. 230l/H
Stirring machine
Speed: 1000min-1 (constant)
Power consumption: 180W.



Laboratory Instrument India

