



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Suspension Production Unit

**Product Code :**  
LBNY-0005-12700034



**Description :**

Eccentric Screw Pump, With Pressure Cut-Out Switch, Dry-Running Protection And Adjustable Speed, To Deliver the suspension  
Supply Unit To Produce And Deliver Suspensions For Experimental Filtration Units  
Stirred Tank With Lid And Stirring Machine To Prepare A Suspension

**Technical Specification :**

Tank: 200L, stainless steel  
Pump  
Manometer measuring range: 0...10bar  
Max. Head: 50m  
Max. Flow Rate: Approx. 230l/H  
Stirring machine  
Speed: 1000min<sup>-1</sup> (constant)  
Power consumption: 180W.



**Laboratory Instrument India**

