



Email : sales@laboratoryinstrumentindia.com

Product Name :
Nutsche Pressure Filter

Product Code :
LBNY-0005-12700037



Description :

Nutsche Pressure Filter For Discontinuous Cake Filtration
Bottom Flange With Recessed Sieve Base And PP Filter Cloth
Enclosed 3-Part Vessel With 2 Flanges And 2 Bumped Bases
Bottom Section Of Vessel To Collect Filtrate
Maintenance And Pressure Control Unit To Adjust Positive Pressure In Centre And Top Section
Centre Section To Form Filter Cake
Top Section Removable To Remove Filter Cake
2 Sight Glasses To Observe Level In Bottom Section
2 Manometers To Indicate Pressure Upstream And Downstream Of Filter
Production and transport of suspension with suspension production unit.

Technical Specification :

Vessel
Material: stainless steel
Capacity: approx. 75l
Inside diameter: approx. 300mm
Permissible temperature: -10...100°C
Permissible pressure: -1...10bar
Measuring ranges
1x maintenance and pr
2x manometers (D=160mm): 0...4bar.



Laboratory Instrument India