



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Nutsche Pressure Filter

**Product Code :**  
LBNY-0005-12700037



#### **Description :**

Nutsche Pressure Filter For Discontinuous Cake Filtration  
Bottom Flange With Recessed Sieve Base And PP Filter Cloth  
Enclosed 3-Part Vessel With 2 Flanges And 2 Bumped Bases  
Bottom Section Of Vessel To Collect Filtrate  
Maintenance And Pressure Control Unit To Adjust Positive Pressure In Centre And Top Section  
Centre Section To Form Filter Cake  
Top Section Removable To Remove Filter Cake  
2 Sight Glasses To Observe Level In Bottom Section  
2 Manometers To Indicate Pressure Upstream And Downstream Of Filter  
Production and transport of suspension with suspension production unit.

#### **Technical Specification :**

Vessel  
Material: stainless steel  
Capacity: approx. 75l  
Inside diameter: approx. 300mm  
Permissible temperature: -10...100°C  
Permissible pressure: -1...10bar  
Measuring ranges  
1x maintenance and pr  
2x manometers (D=160mm): 0...4bar.



**Laboratory Instrument India**