



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Stirring

**Product Code :**  
LBNY-0005-12700038



**Description :**

Transparent Stirred Tank With 4 Removable Flow Impeders  
Investigation Of Mixing Processes During Stirring  
Speed-Controlled Stirring Machine With Digital Torque Indicator  
9 Interchangeable Stirrers:  
Portable device for measuring conductivity and temperature  
Axial-, Radial-, Tangential-Conveying  
Removable Coiled Tube For Cooling Or Heating With External Water Supply.

**Technical Specification :**

Stirring



**Laboratory Instrument India**