



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Fundamentals Of Sedimentation

**Product Code :**  
LBNY-0005-12700039



**Description :**

5 transparent tanks with scale for comparison of the settling velocities of solids in various suspensions  
Experiments in the fundamentals of sedimentation  
3 measuring cups for preparation of suspensions  
Pycnometer to determine the density of the liquids and solids  
Tanks removable for filling, mixing and cleaning  
Recommended accessories: balance, coagulant  
Tanks backlit by fluorescent tubes to aid observation  
Stopwatch to record the sedimentation time.

**Technical Specification :**

Tanks  
Material: PMMA  
Length: 1000mm  
Scale division: 1mm  
Inside diameter: 42mm  
Measuring cups  
Capacity: 2000ml  
Scale division: 50ml  
Fluorescent Tubes  
Power: 6x 18W  
Stopwatch  
Resolution: 1/100sec.



**Laboratory Instrument India**