




Email : sales@laboratoryinstrumentindia.com

Product Name :
Pneumatic Transport

Product Code :
LBNY-0005-12700042



Description :

4 interchangeable injectors to disperse the feed material into the air flow
Vertical tube made of glass
Pneumatic pressure-lifting of solids in a vertical tube
Feed of solid into air flow via vibrating trough with adjustable throw
Measuring points for pressure loss, temperature and flow velocity
Precision pressure regulator to adjust input pressure and flow rate
Collector and feed tanks interconnected by tube with plug valve
Collector and feed tanks made of transparent material (pmma).

Technical Specification :

Vertical tube
Diameter: 50mm
Height: 2m
Tanks
Collector: 40L
Feed: 20L.

Laboratory Instrument India

