



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Rolling Agglomeration

**Product Code :**  
LBNY-0005-12700048



#### **Description :**

Air Pressure Adjustment By Pressure Reducing Valve  
Tanks for solid, granulating liquid and agglomerates  
Rolling Agglomeration With A Dish Granulator  
Positions Of Solid And Liquid Feed Adjustable  
Metering Device To Adjust The Mass Flow Of Solid Feed Material  
Two-Component Nozzle To Atomise The Granulating Liquid With Compressed Air  
Dish Granulator With Adjustable Rotation Speed And Angle Of Inclination  
Peristaltic Pump To Adjust The Flow Rate Of Liquid.

#### **Technical Specification :**

Dish granulator  
Material: stainless steel  
Diameter: Approx. 400mm  
Rim Height: Approx. 100mm  
Pump  
Max. flow rate: approx. 428mL/min  
Dish drive motor  
Speed: 20...400min<sup>-1</sup>  
Power Consumption: Approx. 750w  
Tanks  
Agglomerates: 10l  
Solids: 40L

---

Solids silo: approx. 10l  
Granulating liquid: 5l.



**Laboratory Instrument India**