



Email : sales@laboratoryinstrumentindia.com

Product Name :
Rolling Agglomeration

Product Code :
LBNY-0005-12700048



Description :

Air Pressure Adjustment By Pressure Reducing Valve
Tanks for solid, granulating liquid and agglomerates
Rolling Agglomeration With A Dish Granulator
Positions Of Solid And Liquid Feed Adjustable
Metering Device To Adjust The Mass Flow Of Solid Feed Material
Two-Component Nozzle To Atomise The Granulating Liquid With Compressed Air
Dish Granulator With Adjustable Rotation Speed And Angle Of Inclination
Peristaltic Pump To Adjust The Flow Rate Of Liquid.

Technical Specification :

Dish granulator
Material: stainless steel
Diameter: Approx. 400mm
Rim Height: Approx. 100mm
Pump
Max. flow rate: approx. 428mL/min
Dish drive motor
Speed: 20...400min⁻¹
Power Consumption: Approx. 750w
Tanks
Agglomerates: 10l
Solids: 40L

Solids silo: approx. 10l
Granulating liquid: 5l.



Laboratory Instrument India