



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Osborne Reynolds Experiment

**Product Code :**  
LBNY-0005-12700049



**Description :**

Water tank with glass beads to stabilise the flow  
Flow rate in the pipe section can be adjusted via a valve  
Visualisation of laminar and turbulent flow in the Osborne Reynolds experiment  
Vertical glass pipe section  
Flow rate determined by base module  
Water supply using base module or via laboratory supply  
Water as flowing medium and ink as contrast medium.

**Technical Specification :**

Tank for ink  
Capacity: approx. 250mL  
Water tank  
Capacity: 2200mL  
Pipe section  
Inside diameter: 10mm  
length: 675mm.

---

# Laboratory Instrument India

