



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Assembly Exercise Piston Compressor

**Product Code :**  
LBNY-0005-12700051



**Description :**

Single-stage compressor, air-cooled, with fan flywheel, intake filter and pipe unions  
Assembly exercise for engineering training  
Compressor consisting of piston and cylinder, housing, driving gear, cylinder cover with valves  
Complete, disassembled piston compressor with small parts set and 6 assembly jigs in a storage case.

**Technical Specification :**

Air-cooled single-cylinder piston compressor  
Intake Capacity: 115l/Min  
Displacement Volume: 63cm<sup>3</sup>  
Speed: 1850min<sup>-1</sup>  
Cylinder Bore: 50mm  
Stroke: 32mm  
Max. Pressure: 10bar  
Dimensions, assembled, LxWxH: 223x256x314mm  
Drive Power Output: 700w.



**Laboratory Instrument India**

