



Email : sales@laboratoryinstrumentindia.com

Product Name :
Trainer For Various Heat Exchangers

Product Code :
LBNY-0005-12700055



Description :

Electromagnetic flow meter
Digital displays for temperature, pressure differences and flow rate
Hot and cold water benches available for independent operation
Examination and comparison of various heat exchanger types
Water/steam heat exchanger unit and electrical steam generator available for further experiments
Operating mode (parallel flow or counter flow) selectable via valves
Flow rates adjustable via valves
Five different types of heat exchangers included in the scope of delivery
Finned heat exchanger with fan.

Technical Specification :

Trainer For Various Heat Exchangers



Laboratory Instrument India