



Email : sales@laboratoryinstrumentindia.com

Product Name :
Fundamentals Of Humidity Measurement

Product Code :
LBNY-0005-12700056



Description :

Humidification via ultrasonic atomiser
Dehumidification via Peltier cooling element
Different measuring methods for measuring humidity
Climatic chamber with adjustable humidity and transparent door
Fan for air recirculation
2 electronic instruments: capacitive sensor, hygrometer with synthetic fibre and combined temperature sensor
2 mechanical instruments: psychrometer, hair hygrometer.

Technical Specification :

Humidifier
Low water cut-off
Power consumption: 21,6W
Ultrasonic atomiser

Dehumidifier
Peltier element
Cooling surface: 1600mm²
Cooling capacity: 56,6W (50°C ambient temperature).

Laboratory Instrument India

