



Email : sales@laboratoryinstrumentindia.com

Product Name :
Rotary Vane Vacuum Pump

Product Code :
LBNY-0005-12700057



Description :

Oil separator with tank on the delivery side of the rotary vane pump
Vacuum generation with a rotary vane pump
Manometer to display the pressure in the vessel
2 silencers to reduce flow noise
1 inlet line with needle valve, 1 inlet line with needle valve and flow meter on the pressure vessel for Simulation loads and leaks.
Vacuum generation in a steel pressure vessel.

Technical Specification :

Rotary vane pump
Final vacuum: 0,02mbar
Motor power output: 370W
Volume of pressure vessel: 20L
Max. flow rate: 5m³/h
Oil separator
Capacity of tank: 150mL
Degree of separation: 99,98%.

Laboratory Instrument India

