



Email : sales@laboratoryinstrumentindia.com

Product Name :
Assembly Exercise Piston Compressor Functional Test

Product Code :
LBNY-0005-12700059



Description :

Pressure vessel with adjustable pressure switch and manometer
Switch box with controls and ammeter to indicate current consumption
Unit for functionality testing of the piston compressor from assembly exercises
Safety devices: pressure switch, safety valve, protective screen for belt drive, emergency-off switch
Driven by electric motor and belt drive
Single-phase electric motor on adjustable carriage.

Technical Specification :

Pressure vessel
Max. pressure: 10bar
Capacity: 10L
Single-phase motor
Speed: 1405min⁻¹
Power output: 250W
Measuring ranges
Ammeter: 0...4A, class 2,5
Manometer: 0...16bar.

Laboratory Instrument India

