



Email : sales@laboratoryinstrumentindia.com

Product Name :
Assembly Maintenance Exercise Diaphragm Pump

Product Code :
LBNY-0005-12700060



Description :

Flow setting by manual stroke length adjustment (including during operation)
Manual drive with crank instead of a drive motor
Learning concept for maintenance and repair exercises on a single-diaphragm pump
The kit forms part of the GUNT assembly, maintenance and repair practice line
Diaphragm and push rod directly linked
Pump parts and tools housed in a tool box.

Technical Specification :

Assembly Maintenance Exercise Diaphragm Pump



Laboratory Instrument India