



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Assembly Maintenance Exercise Diaphragm Pump

**Product Code :**  
LBNY-0005-12700060



**Description :**

Flow setting by manual stroke length adjustment (including during operation)  
Manual drive with crank instead of a drive motor  
Learning concept for maintenance and repair exercises on a single-diaphragm pump  
The kit forms part of the GUNT assembly, maintenance and repair practice line  
Diaphragm and push rod directly linked  
Pump parts and tools housed in a tool box.

**Technical Specification :**

Assembly Maintenance Exercise Diaphragm Pump



**Laboratory Instrument India**