



Email : sales@laboratoryinstrumentindia.com

Product Name :
Multi Channel Measuring Amplifier

Product Code :
LBNY-0005-12700061



Description :

Amplifier
Strain gauge connection via 6-pin single ports or 68-pin multiports
Extension of measuring points possible
Multi-channel measuring amplifier for processing of strain gauge signals
Processing of measured values on a PC
Strain gauge connection in half or full bridge configuration
Connection to the PC via USB.

Technical Specification :

Amplifier
Max. 32 input channels
Number of input channels: 16.



Laboratory Instrument India