



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Centripetal Force Pendulum

**Product Code :**  
LBNY-0004-780025



**Description :**

Centripetal Force Pendulum

**Technical Specification :**

Centripetal Force pendulum:

The Centripetal Force Pendulum attaches to a Force Sensor and allows students to directly measure the forces involved in circular motion. By attaching the Force Sensor/pendulum combination to the Rotary Motion Sensor, the relationship between force, mass, and velocity in a circular path can be investigated.

Measures the force vs the angular velocity

Include pendulum rod with threaded connector

Sliding mass and a cord clip



**Laboratory Instrument India**