



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Centrifugal Pump Standard Design

**Product Code :**  
LBNY-0005-12700065



**Description :**

Powered by Universal Drive and Brake Unit  
Pressure sensors at the inlet and outlet of the pump  
Examination of a standard centrifugal pump  
Pressure display on the display unit  
Operation with Supply Unit for Water Pumps

**Technical Specification :**

Standard centrifugal pump  
Max. head: 22m  
Nominal speed: approx. 2900min<sup>-1</sup>  
Max. flow rate: 24m<sup>3</sup>/h.



**Laboratory Instrument India**