



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Two Stage Piston Compressor

**Product Code :**  
LBNY-0005-12700067



**Description :**

Pressure vessel after the second stage with safety valve, blow-off valve and silencer as well as an Additional manometer and a pressure switch  
Intake vessel with nozzle to measure the intake volumetric flow rate, pressure sensor and additional manometer  
Recording the characteristic of a two-stage compressor  
Digital displays for temperatures, pressures, differential pressures and electric power output  
Piston compressor with 2 cylinders in V-arrangement  
Pressure vessel after the first stage as intercooler  
Sensors for pressures, temperatures and electric power output.

**Technical Specification :**

Two Stage Piston Compressor



**Laboratory Instrument India**