



Email : sales@laboratoryinstrumentindia.com

Product Name :
Multi Head Diaphragm Pump

Product Code :
LBNY-0005-12700068



Description :

Common stroke rate adjustable for all 3 heads
Stroke individually adjustable for each head
Examination of a multi-head diaphragm pump
Diaphragm pump with 3 heads
3 scaled tanks to determine the flow rate
1 feed tank
6 manometers to measure the pressure on the intake and delivery side
3 overflow valves to set the counter pressure.

Technical Specification :

3 overflow valves
Adjustable pressure range: each 2...6bar
Tanks
Tanks with scale: each 5L
Feed tank: 60L
Pump
Max. stroke: each approx. 6mm
Flow rate: 3x 0...18L/h
Max. head: each approx. 100m
Max. suction head: each approx. 3m.



Laboratory Instrument India