



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Multi Head Diaphragm Pump

**Product Code :**  
LBNY-0005-12700068



**Description :**

Common stroke rate adjustable for all 3 heads  
Stroke individually adjustable for each head  
Examination of a multi-head diaphragm pump  
Diaphragm pump with 3 heads  
3 scaled tanks to determine the flow rate  
1 feed tank  
6 manometers to measure the pressure on the intake and delivery side  
3 overflow valves to set the counter pressure.

**Technical Specification :**

3 overflow valves  
Adjustable pressure range: each 2...6bar  
Tanks  
Tanks with scale: each 5L  
Feed tank: 60L  
Pump  
Max. stroke: each approx. 6mm  
Flow rate: 3x 0...18L/h  
Max. head: each approx. 100m  
Max. suction head: each approx. 3m.



**Laboratory Instrument India**