



Email : sales@laboratoryinstrumentindia.com

Product Name :
Controlled System Module Level

Product Code :
LBNY-0005-12700071



Description :

Transparent glass tank with overflow, drain valve and manometer
Glass tank can be hermetically sealed
Level control loop setup
Adjustable safety valve
Fitting for attaching pressure sensor or level sensor
Ready-to-install compact panel assembly
Plexiglass protective sleeve
PVC tank base and cover
Cascade control option in conjunction with RT 450.02
Water as the working medium, supply via base module RT 450
Plexiglass protective sleeve
PVC tank base and cover.

Technical Specification :

Tank
Capacity: max. 7L
Manometer: 02,5bar
Safety valve: 1...4bar adjustable
Filling level: max. 475mm
Operating pressure: max. 2bar.



Laboratory Instrument India