



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Projectile Launcher

**Product Code :**  
LBNY-0004-780026



**Description :**

Projectile Launcher

**Technical Specification :**

Mechanical Projectile accessory 01:

Projectile launcher:

Adjustable incline and power

Projectile range: up to 5 meters

Include a stand to fix the launcher with a pack of large C-Clamp to fixed on the table

Provide necessary accessories:

Angle indicator

? 3 projectiles (plastic balls or equivalent)

1 plunger and

? 2 safety goggles

Item also include:

Photogate mounting bracket

C-clamp to fix the launcher on the table, at least qty 6

Set of two photogate heads

Large base of 4 KG to keep it steady

Multi clamp

Rod of 1 meter length

Rod of 40 cm length



**Laboratory Instrument India**