



Email : sales@laboratoryinstrumentindia.com

Product Name :
Pressure Control Trainer

Product Code :
LBNY-0005-12700073



Description :

2 pressure tanks with pressure relief valve and manometer for direct observation of the tank pressure
Valves permit investigation of a 1st order controlled system (1 tank) or 2nd order controlled system (2 in-line tanks)
Trainer for control engineering experiments
Generation of disturbance variables by drain valve
Pneumatically operated control valve with electro-pneumatic positioned digital controller, Parameterisable as a P, PI or PID controller
Pressure control process, equipped with standard industrial components
Pressure measurement by pressure sensor
Process variables X and Y accessible as analogue signals via lab jacks.
2-channel line recorder.

Technical Specification :

Pressure Control Trainer



Laboratory Instrument India

