



Email : sales@laboratoryinstrumentindia.com

Product Name :
Process Automation Training System Base Module

Product Code :
LBNY-0005-12700074



Description :

Practice-oriented terminal blocks for wiring of the instrumentation and control signals
2x maintenance unit and a pressure regulator for compressed air
Frame for electrical modules
Water circuit with centrifugal pump and tank
Base module for a flexible process automation training system
Large frame with aluminium rails
Switch cabinet with 24V power supply
Installation and connecting material including tools for connecting accessory modules.

Technical Specification :

24VDC supply for components
Power consumption: 370W.
Stainless steel storage tank: 75L
Centrifugal pump
Max. head: 20m
Max. flow rate: 5,4m³/h.

Laboratory Instrument India

