



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Process Automation Training System Base Module

**Product Code :**  
LBNY-0005-12700074



**Description :**

Practice-oriented terminal blocks for wiring of the instrumentation and control signals  
2x maintenance unit and a pressure regulator for compressed air  
Frame for electrical modules  
Water circuit with centrifugal pump and tank  
Base module for a flexible process automation training system  
Large frame with aluminium rails  
Switch cabinet with 24V power supply  
Installation and connecting material including tools for connecting accessory modules.

**Technical Specification :**

24VDC supply for components  
Power consumption: 370W.  
Stainless steel storage tank: 75L  
Centrifugal pump  
Max. head: 20m  
Max. flow rate: 5,4m<sup>3</sup>/h.

---

## Laboratory Instrument India

