



Email : sales@laboratoryinstrumentindia.com

Product Name :
Controlled System Module Flow

Product Code :
LBNY-0005-12700077



Description :

Flow control loop setup
Stop valve in the inlet section
Water as the working medium, supply via base module
Ready-to-install compact panel assembly
Optional cascade control in conjunction
Rotameter with scale.

Technical Specification :

Rotameter: 0...2,5m³/h
Dimensions And Weight
LxWxH: 250x180x700mm
Weight: approx. 10kg.



Laboratory Instrument India