



Email : sales@laboratoryinstrumentindia.com

Product Name :
PLC Module

Product Code :
LBNY-0005-12700080



Description :

PLC with 2 integrated set point encoders
Module for basic exercises on a programmable logic controller (PLC)
Example program supplied
Integrated patch board for creating circuits with input and output elements
Self-contained PLC module, usable as part of a complex system.

Technical Specification :

Memory type: PLC back-up battery for 32kByte RAM and clock
Rated voltage: 24VDC.
Connections
2 analogue inputs
1 analogue output
16 digital inputs
16 digital outputs.



Laboratory Instrument India

