



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
PLC Module

**Product Code :**  
LBNY-0005-12700080



#### **Description :**

PLC with 2 integrated set point encoders  
Module for basic exercises on a programmable logic controller (PLC)  
Example program supplied  
Integrated patch board for creating circuits with input and output elements  
Self-contained PLC module, usable as part of a complex system.

#### **Technical Specification :**

Memory type: PLC back-up battery for 32kByte RAM and clock  
Rated voltage: 24VDC.  
Connections  
2 analogue inputs  
1 analogue output  
16 digital inputs  
16 digital outputs.



**Laboratory Instrument India**

