



Email : sales@laboratoryinstrumentindia.com

Product Name :
Controlled System Module Pressure

Product Code :
LBNY-0005-12700085



Description :

Construction of a pressure control loop
2 pressure tanks with safety valves
Direct pressure indication by 2 manometers
Ready-to-install compact panel assembly
Compressed air as working medium, supply via base module
Valve permits series configuration of both pressure tanks
2 valves: 1x drain, 1x vent
Copper piping
Series configuration of pressure tanks: investigation of coupled controlled system response compared Single-tank operation.

Technical Specification :

Pressure tank
Capacity: 3L
Manometer: 0...10bar
Operating pressure: max. 6bar
Safety valve: adjustable to max. 10bar.

Laboratory Instrument India

