



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Controlled System Module Pressure

**Product Code :**  
LBNY-0005-12700085



**Description :**

Construction of a pressure control loop  
2 pressure tanks with safety valves  
Direct pressure indication by 2 manometers  
Ready-to-install compact panel assembly  
Compressed air as working medium, supply via base module  
Valve permits series configuration of both pressure tanks  
2 valves: 1x drain, 1x vent  
Copper piping  
Series configuration of pressure tanks: investigation of coupled controlled system response compared Single-tank operation.

**Technical Specification :**

Pressure tank  
Capacity: 3L  
Manometer: 0...10bar  
Operating pressure: max. 6bar  
Safety valve: adjustable to max. 10bar.

---

# Laboratory Instrument India

