



Email : sales@laboratoryinstrumentindia.com

Product Name :
Time of Flight Accessory

Product Code :
LBNY-0004-780028



Description :

Time of Flight Accessory

Technical Specification :

Mechanical projectile accessory 03:

Time of flight accessory.

The Time-of-Flight Accessory is designed primarily for freefall or projectile motion experiments. When an object hits the plate, a signal is sent to the interface.

This accessory is used with the projectile launcher to calculate the time of flight of the ball launched from the projectile.

Include instruction manual and experiment guide.

Includes: extension cable.



Laboratory Instrument India