



Email : sales@laboratoryinstrumentindia.com

Product Name :
Networking Of Industrial Controllers

Product Code :
LBNY-0005-12700093



Description :

All process variables accessible as analogue signals at lab jacks
Experimental unit for networking of industrial controllers
2 signal generators
Profibus DP interface card for PC
2 digital controllers, configurable as P, PI or PID controllers, with field bus interface.

Technical Specification :

Controller
Integral action time T_n : 0...3.600s
Signal generator: 0...10V
Derivative time T_v : 0...1.200s
Configurable as P, PI or PID controller
Proportional gain X_p : 0...999,9%
Process variables as analogue signals: 0...10V
Connection of external instruments (e.g. oscilloscope, line recorder) via lab jacks.



Laboratory Instrument India

