



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Calibration Station

**Product Code :**  
LBNY-0005-12700094



**Description :**

Sending and measuring pneumatic and electrical signals  
3 pressure regulators with manometers  
Adjustable power supply unit as DC voltage source  
Investigation of transmission behaviour and calibration of electronic and pneumatic control loop components  
3 digital pressure gauges with differing measuring ranges  
2 adjustable potentiometers  
2 switchable AC voltage sources with differing voltages  
3-channel line recorder with freely selectable measuring ranges  
Separate high-grade universal calibrator for voltage, current, temperature and resistance  
Freely configurable industrial controller  
Digital voltmeter  
Digital ammeter.

**Technical Specification :**

Calibration Station



**Laboratory Instrument India**

