



Email : sales@laboratoryinstrumentindia.com

Product Name :
Calibration Station

Product Code :
LBNY-0005-12700094



Description :

Sending and measuring pneumatic and electrical signals
3 pressure regulators with manometers
Adjustable power supply unit as DC voltage source
Investigation of transmission behaviour and calibration of electronic and pneumatic control loop components
3 digital pressure gauges with differing measuring ranges
2 adjustable potentiometers
2 switchable AC voltage sources with differing voltages
3-channel line recorder with freely selectable measuring ranges
Separate high-grade universal calibrator for voltage, current, temperature and resistance
Freely configurable industrial controller
Digital voltmeter
Digital ammeter.

Technical Specification :

Calibration Station



Laboratory Instrument India

