



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Calibration Trainer

**Product Code :**  
LBNY-0005-12700095



**Description :**

2 DC voltage sources to supply control loop components with auxiliary power  
2 precision manometers  
2 pressure regulators with manometers  
Height-adjustable tank and scale to set low pressures  
2 controllers to send and measure current signals  
Investigation of transmission behaviour and calibration of control loop components and measuring instruments  
Sending and measuring pneumatic and electrical signals  
Control loop components and meters available as accessories.

**Technical Specification :**

2 pressure regulators  
0...6bar  
0...1,6bar  
Pressure range, height-adjustable tank  
0...1000mm head  
1 output each: 4...20mA  
1 input each: 4...20mA  
2 controllers  
2 DC voltage sources  
each 24VDC  
2 precision manometers  
0...6bar, D=160mm

---

0...1,6bar, D=160mm.



**Laboratory Instrument India**